

Quanta System

Chrome





THE ULTIMATE
EXPANDABLE
LASE STATION

Chrome[®]

L A S E - S T A T I O N

The ultimate **expandable Lase Station** based on a **reinvented Q-Switched technology** to maximize the operator possibilities and quality of outcomes.



DRIVE THE CHANGE

Combining strong innovation with maximum flexibility increases consumers' involvement in treatment selection.

More customer satisfaction and retention.



EXPAND YOUR POSSIBILITIES

100% Italian design with 35 years of heritage in the development of lasers for people care.



SELECT QUANTA SKIN APPS[®]

Chrome Lase Station empowers you to perform more than 50 kind of treatments thanks to additional modules and exclusive Skin Apps[®].




INCREASE YOUR BUSINESS

The Lase Station that grows together with your business. Upgrade anytime with the option you need to enhance your practice.

Drive Your — Change

Finally you will have in your hands a never-seen-before platform empowering you to offer **more than 50** different kind of **treatments** and the exclusive **Quanta Skin Apps®** designed to explain in a contemporary way new aesthetic solutions to your consumers.





Fractional HC

Fractional DF

Nd:YAG
1064 nm P

TWAIN
IPL

ROSSO
694 nm QS

Nd:YAG
532 nm
1064 nm
QS - OP

TWAIN
Er:YAG 2940

NaturaPeel



Expand Your — Possibilities

Chrome Lase Station empowers you to perform **a wide range of treatments** thanks to additional **modules**, 100% Italian design with 35 years of **heritage** in the **development of lasers** for people care.

Expand your possibilities

OPTIZOOM

Ultimate optics for the highest treatment precision

Sources available: 1064 Pulsed (optional) | 1064 Q-Switched | 532 Q-Switched



Dermatological
Vascular Lesions



Benign Pigmented
Lesions



Skin
Rejuvenation



Multicolour Tattoo
Removal



Onychomycosis

ROSSO

Ultra durable active 694 nm Q-Switched laser handpiece



Benign Pigmented
Lesions



Tattoo
Removal




Expand your possibilities

NATURAPEEL

Fixed handpiece for cream-assisted laser peeling



NaturaPeel®



Fractional HC

High coverage perfection for fractional treatment of pigmentations

- Skin Rejuvenation
- Benign Pigmented Lesions & Toning



Twain IPL

Mix of Power and versatility

- Benign Vascular Lesions • Benign Pigmented Lesions • Skin Rejuvenation
- Hair Removal



Twain 2940

The ablative power in your hands

- Scars • Ablative Dermatological Procedures • Skin Resurfacing

Fractional DF

Deep fractional multi-level solution for skin resurfacing


- Scars • Skin resurfacing • Skin Rejuvenation • Resistant Tattoo Removal



Select Your — Skin Apps®

Your consumers will receive, from now on, an easy-to-understand offer. This will help you **increase** their personal **involvement** and higher **retention**.







NaturaPeel®

FAST AND GENTLE LASER PEEL

The laser pass on the natural cream, based on bamboo carbon extracts, removes the superficial impurities.



NATURAL LASER PEEL CREAM

The Natural cream makes the laser peel treatment unique preserving the benefits from the original cream with a more eco-friendly and sustainable procedure.



TREAT-AND-GO

The treatment is free of down-time with immediate effect of smoothness and reduced facial pore sizes.



DÉCOLLETÉ

REDUCTION OF REDNESS

For a décolleté skin free from undesired skin vascularization caused by sun exposure.



tone uniformity

For a more homogenized décolleté skin free from undesired photo-aged dark spots and age signs.



SKIN TIGHTENING EFFECT

To promote neocollagenesis helping the skin being more tighten with wrinkles reduction.



NATURAPEEL

Courtesy by Quanta System



DÉCOLLETÉ

Courtesy by Quanta System



Refresh

TONE UNIFORMITY

For a more homogenized skin free from undesired pigmentations and age signs.



SKIN BRIGHTNESS

For a real refresh from the damages of working fatigue and unprotected sunbathing.




SKIN REFILL

To gently recharge consumers' skin giving it the right energy to tackle the ageing process.



REFRESH

Courtesy by Nicola Zerbinati MD



Refresh Bronze

for darker skin types

THE REFRESH ADVANTAGES

The same REFRESH Skin App clinical outcomes.



DARKER SKINS

This option has been developed to provide dark skin consumers with the right attention they deserve.



REFRESH BRONZE

Courtesy by Nicola Zerbinati MD



Define

TEXTURE RE-DEFINITION

To remove the skin irregularities and imperfections caused by aging processes.



LASER PEELING

To perform a gentle and non-invasive skin peeling to promote a faster epidermal turn-over for a younger appearance.



SKIN RESURFACING

To promote the old skin replacement with dermal remodelling by deep neocollagenesis cascade.



Revivo

SKIN RESURFACING & PEELING

To combine the effects of skin resurfacing and peeling for the rebirth of consumers' skin.



TIGHTENING EFFECT

To restore the original collagen configuration for excellent tissue definition and firming.



COUPEROSE REDUCTION

To efficiently remove undesired vascularization restoring the native condition of the skin.











DEFINE

Courtesy by Nicola Zerbinati MD

REVIVO

Courtesy by Nicola Zerbinati MD

...and much more!

		CLASS OF TREATMENT			
					
		Benign Pigmented Lesions	Benign Cutaneous Lesions	Dermatological Vascular Lesions	Hair Removal
HI-TECH MODULES	Basic	●	●		
	Pulsed	●	● ²	●	●
	NaturaPeel				
	Fractional DF	●	●		
	Fractional HC	●			
	Rosso	●			
	Twain IPL	●	● ³	●	●
	Twain 2940		●		
HI-TECH MODULES					
		Tattoo Removal	Skin Resurfacing	Skin Rejuvenation	Skin Peeling ¹
	Basic	●		●	
	Pulsed			●	
	NaturaPeel			●	●
	Fractional DF	●	●	●	
	Fractional HC			●	
	Rosso	●			
Twain IPL			●		
Twain 2940		●	●		

BASIC: 532/1064 nm Q-Switched Nd:YAG laser

PULSED: 1064 nm pulsed mode based on a Nd:YAG laser source

NATURAPEEL: dedicated handpiece to laser skin peeling

FRACTIONAL DF: deep fractional handpiece

FRACTIONAL HC: high coverage fractional handpiece

ROSSO: active 694 nm Q-Switched laser handpiece (not dye)

TWAIN IPL: powerful IPL handpiece with up to 6 kind of range of wavelengths

TWAIN 2940: full spot and fractional Erbium handpiece


¹ Skin peeling is intended as the result of laser skin rejuvenation of the most superficial epidermal lesions

² Treatment of Onychomycosis and warts

³ Treatment of acne vulgaris and inflammatory acne


Increase Your — Business

The Lase Station that **grows together** with your business. **Expand** anytime with the modules you need to **enhance your practice**.



SCALABLE STATION

The Chrome Lase Station® can be expanded just with the addition of the laser modules that are needed.



MORE THAN 50 INDICATIONS

A wide range of indications is available for a great number of treatments and therapies.



POWERFUL ACCESSORIES

Multiple powerful accessories and laser sources can be added for a complete setup of configuration.



RESULT ORIENTED

The Lase Station is totally dedicated to get the maximum clinical results and consumers' satisfaction.



GREAT APPEAL

Unique and exclusive Skin Apps® to efficiently communicate new treatments possibilities to take care of your consumers' skin.

Technical Specifications

SOURCES	QS/OP Nd:YAG	QS/OP Nd:YAG	P Nd:YAG (optional)
Wavelength (nm)	1064	532	1064
Pulse Width	6ns in QS 6ns + 6ns in OP	6ns in QS 6ns + 6ns in OP	0.3 ÷ 50 ms
Burst Mode ¹	N/A		Up to 3
Spot Size (mm) ²	2 ÷ 8 mm round NaturaPeel 7 mm round 8 mm DF fractional round 9 mm HC fractional round 2x2, 3x3, 4x4, 5x5, 7x7 mm ² square		
Max Fluence (J/cm ²)	32 (QS) 48 (OP)	16 (QS) 19 (OP)	500 (SW Limited)
Repetition Rate (Hz)	Up to 20		Up to 10 (step 0.5)
Legend	QS Q-Switched OP Opti-Pulse P Pulsed		

¹ Not available for rep. rates higher than 4 Hz.

² The availability of the mentioned handpiece spot sizes is subject to the marketing authorization in some Country, contact Sales Repr.

OPTION	Twain 2940	Twain IPL						Rosso
Wavelength (nm)	2940	650-1200	625-1200	590-1200	570-1200	550-1200	400-1200	694
Pulse width	0.3 ÷ 1ms	3 - 40ms (pulses train duration up to 240ms)						4 ns
Max Fluence (J/cm ²)	up to 95	Fluence up to 25 (both LARGE and SMALL waveguide)						Up to 5
Repetition rate (Hz)	up to 6	up to 3 Hz						up to 10
Spot size (mm)	2, 4, 9 round; 9 fractional round	48 x 13 (LARGE) 25 x 13 (SMALL)						2, 3, 4 round

Quanta System

This brochure is not intended for the U.S. market.
Certain Intended Uses/Configurations/Models/Accessories are not cleared for U.S.



*Nayel Medical is the exclusive distributor for
Quanta System Aesthetic Lasers in UAE*

DANGER - VISIBLE AND INVISIBLE LASER RADIATION

Avoid eye or skin exposure to direct or scattered radiation | Laser product: Class 4 | Aiming beam: Class 3R

Note: National local authorities may put restrictions to the parameters indicated in the table in the previous page, or may limit or remove certain intended uses. Specifications are subject to change without notice.



Quanta System S.p.A. – Via Acquedotto, 109, 21017 Samarate (VA) – Italy
© 2023 All rights reserved – quanta@quantasystem.com | www.quantasystem.com