

Quanta System

UNLEASH THE
◊ SUPREMA ◊
POWER

 NAYEL MEDICAL

*Nayel Medical is the exclusive distributor
for Quanta System Aesthetic Lasers in UAE*

◆ SUPREMA ◆

Suprema is the new laser platform developed by Quanta System, dedicated to vascular treatments, hair removal and rejuvenation with fast and effective results and reduced down-time. This technology is among the highest performing in the world thanks to its 3 wavelengths (532 nm, 755 nm and 1064 nm) and exclusive **Quanta Mixed Technology®**.

◆ Handpieces

Suprema has zoom optics (2: 2-6 mm and 8-20 mm) and fixed optics (14,16,18,20,22 mm) that can be easily adapted to the size of the areas to be treated.

One-click Connection

This system automatically detects spots without external wiring.

◆ Skin Cryo & Skin Cooler

Cooling systems, revised to their most ergonomic form, maximize treatment comfort.

New User Experience

The intuitive graphical user interface GUI and unprecedented treatment database complement the laser's unique user experience.

Twain Handpieces

Twain Er:YAG 2940 nm or Twain IPL handpieces can be connected to further extend the indications for use.

Quanta Mixed Technology®

The technology that allows, simultaneously or sequentially, to mix and emit 2 wavelengths: Alexandrite (755 nm) and Nd:YAG (1064 nm).

Ideal for Vascular and Hair Removal Treatments

Thanks to the different wavelengths, you can act deep into the blood vessels or follicle, with optimal results on face and body, for all skin types.



DANGER - VISIBLE AND INVISIBLE LASER RADIATION Avoid eye or skin exposure to direct or scattered radiation | Laser product: Class 4