



# LITHOBOX<sup>®</sup> & LITHOBOX ZERO<sup>®</sup> INTRACORPOREAL LITHOTRIPTERS

TREATMENT  
PNL  
URS  
TREATMENT



ALL PNEUMATIC, ULTRASONIC PROBES AND HANDPIECES  
ARE COMPATIBLE WITH EMS BRAND DEVICES





# KIDNEY STONE TREATMENT IS MUCH MORE EASIER WITH **LITHOBOX ZERO**





- > All in one
- > Effective
- > Secure
- > Ergonomic Design
- > Portable

*You will see the difference when you use it.*



Touch Screen Interface  
Lithobox Zero



TREATMENT  
PNL  
URS  
TREATMENT





LITHOBOX ZERO is an Intracorporeal Ultrasonic Pneumatic Combined Kidney Stone Crushing System. This design provides a special advantage in combination with clinical data. It is based on the principle of re-absorption with the disintegration of the stone and its aspiration feature.

Lithobox ZERO is designed for combined use. The system enables the operation to be completed much easier and faster by means of pneumatic probes while crushing stones ultrasonically.

- > An easier and faster operation with combined use
- > Longer life different structures ultrasonic and pneumatic probes
- > Ergonomic choice that does not tire the hand
- > Strong suction feature





# LITHOBOX




## PROVIDES FAST, EFFECTIVE AND SAFE RESULTS



LED screen  
With Built-in Compressor  
Lithobox

- > Effective and successful wide usage area
- > Easy and convenient use
- > Superior quality and economical price
- > 5 times stronger closed system and locking pneumatic handle
- > 5 different probe options

*Lithobox offers more than what you expect from a pneumatic lithotripter*





Pneumatic Lithotripter System works on the principle of electro-pneumatic ballistic energy applied by the impact of the projected air by the impact of the projectile on the kidney stone.


Pneumatic lithotripter, is the most effective and powerful method for quickly break the stones have all kinds of density.

- > 100% waterproof closed system
- > It does not require any system for sterilization
- > An easier and faster operation with pneumatic use
- > Longer life different structures pneumatic probes
- > Ergonomic choice that does not tire the hand
- > Extremely stable, steady and seamlessly






**> Handpiece**





**PNEUMATIC HANDPIECE** ●  
Order No.  
IMPL-702







**ULTRASONIC HANDPIECE** ●  
Order No.  
IMPK-601







●  
**PNEUMATIC & ULTRASONIC HANDPIECE**





















<p><b>PNEUMATIC PROBE SET</b></p>  <p>Order No. <b>IMPL-725</b></p> <p><b>SET 0,8/1,0/1,2/1,5/2,0 mm Dia</b></p>	<p><b>PNEUMATIC PROBE</b></p>  <p>Order No. <b>IMPL-715</b></p> <p><b>0,8 mm (2,4 Fr) dia, 640 mm Length</b></p>
<p><b>PNEUMATIC PROBE</b></p>  <p>Order No. <b>IMPL-716</b></p> <p><b>1,0 mm (3,0 Fr) dia, 640 mm Length</b></p>	<p><b>PNEUMATIC PROBE</b></p>  <p>Order No. <b>IMPL-717</b></p> <p><b>1,2 mm (3,6 Fr) dia, 640 mm Length</b></p>
<p><b>PNEUMATIC PROBE</b></p>  <p>Order No. <b>IMPL-718</b></p> <p><b>1,5 mm (4,5 Fr) dia, 640 mm Length</b></p>	<p><b>PNEUMATIC PROBE</b></p>  <p>Order No. <b>IMPL-719</b></p> <p><b>2,0 mm (6,0 Fr) dia, 450 mm Length</b></p>



<p><b>PNEUMATIC PROBE</b></p>  <p>Order No. <b>IMPL-720</b></p> <p><b>3,0 mm (9,0 Fr) dia, 640 mm Length</b></p>	<p><b>ULTRASONIC PROBE SET</b></p>  <p>Order No. <b>IMPK-685</b></p> <p><b>SET (1,5/3,3/3,3/3,8/3,8 mm Dia)</b></p>
<p><b>ULTRASONIC PROBE</b></p>  <p>Order No. <b>IMPL-603</b></p> <p><b>1,5 mm (4,5 Fr) dia, 575 mm Length</b></p>	<p><b>ULTRASONIC PROBE</b></p>  <p>Order No. <b>IMPK-604</b></p> <p><b>3,3 mm (9,9 Fr) dia, 403 mm Length</b></p>
<p><b>ULTRASONIC PROBE</b></p>  <p>Order No. <b>IMPK-605</b></p> <p><b>3,8 mm (11,4 Fr) dia, 403 mm Length</b></p>	<p><b>LITHOBOX ZERO COMBINATION PROBE</b></p>  <p>Order No. <b>IMPK-686</b></p> <p><b>(1,0 mm (3,0 Fr) Dia, 570 mm Length)</b></p>














<p><b>PNEUMATIC HANDPIECE</b></p>  <p>Order No. IMPL-702 LITHOBOX &amp; LITHOBOX ZERO</p>	<p><b>ULTRASONIC HANDPIECE</b></p>  <p>Order No. IMPK-601 LITHOBOX ZERO</p>
 <p>Order No. IMPK-683 LITHOBOX ZERO TROLLEY</p>	 <p>Order No. IMPK-606 FOOT PEDAL</p>
 <p>Order No. IMPK-613 SUCTION JAR (Glass)</p>	 <p>Order No. IMPK-614 SUCTION JAR BASKET</p>

 <p>Order No. IMET-875 POWER CORD</p>	 <p>Order No. IMPK-691 COMPRESSOR HOSE</p>
 <p>Order No. IMPL-716 COMBINATION APPARATUS</p>	 <p>Order No. IMPL-717 PNEUMATIC HANDPIECE ALUMINUM CAP</p>
 <p>Order No. IMPK-655 SUCTION JAR INLET HOSE</p>	 <p>Order No. IMPK-692 SUCTION JAR OUTLET HOSE</p>







 <p>Order No. <b>IMPK-693</b></p> <p><b>SUCTION APPARATUS FOR LITHOBOX</b></p>	 <p>Order No. <b>IMPL-731</b></p> <p><b>PNEUMATIC FOOT PEDAL</b></p>
 <p>Order No. <b>IMPL-718</b></p> <p><b>MINI MEDICAL COMPRESSOR WITH TANK</b></p>	 <p>Order No. <b>IMET-817</b></p> <p><b>LITHOBOX TROLLEY</b></p>
 <p>Order No. <b>IMPL-728</b></p> <p><b>LITHOBOX TROLLEY</b></p>	 <p>Order No. <b>IMET-817-1</b></p> <p><b>LITHOBOX TROLLEY</b></p>


<p><b>PNEUMATIC HANDPIECE</b></p>  <p>Order No. <b>IMPK-690</b></p> <p><b>PNEUMATIC HANDPIECE (SWISS LITHOCLAST COMPATIBLE)</b></p>	<p><b>ULTRASONIC HANDPIECE</b></p>  <p>Order No. <b>IMPK-688</b></p> <p><b>ULTRASONIC HANDPIECE(SWISS LITHOCLAST MASTER COMPATIBLE)</b></p>
 <p>Order No. <b>IMPK-689</b></p> <p><b>EMS COMPATIBLE CONNECTOR FOR PNEUMATIC HANDPIECE</b></p>	 <p>Order No. <b>IMPK-690</b></p> <p><b>EMS COMPATIBLE CONNECTOR FOR ULTRASONIC HANDPIECE</b></p>
 <p>Order No. <b>IMÇ-101</b></p> <p><b>CARRYING BAG</b></p>	

**ULTRASONIC  
PNEUMATIC  
COMBINED**





**LITHOBOX &  
LITHOBOX ZERO**  
INTRACORPOREAL  
LITHOTRIPTERS



*inceler medikal*

**İNCELERMEDİKAL.COM**

İnceler Medikal  
İvedik OSB Mah. 1472. Cad.  
No: 120 Yenimahalle  
Ankara/Türkiye

tel: +90 549 255 3345  
fax: +90 312 255 3347  
sales@incelermedikal.com

